


DEQAS January 2017 - 1,25-dihydroxyvitamin D Method Means * (+/- 1 SD)



	Sample 396	Sample 397	Sample 398	Sample 399	Sample 400
■ DiaSorin Liaison XL (n=90)	100.6	73.4	120.8	91.3	107.9
■ IDS EIA (n=8)	115.8	83.9	135.6	120.6	113.8
■ IDS RIA (n=18)	130.4	103.9	144.7	132.3	122.7
■ IDS iSYS (n=23)	113.9	84.3	120.8	114.3	100.6
■ IDS iSYS New (n=7)	93.3	64.8	88.0	89.7	82.9
■ LC-MS/MS (n=12)	101.9	90.1	127.1	106.3	114.8
■ ALTM (n=165)	105.7	79.2	122.7	101.2	107.4