


## DEQAS October 2017 - 1,25-dihydroxyvitamin D Method Means \* ( $\pm$ 1 SD)



	Sample 411	Sample 412	Sample 413	Sample 414	Sample 415
■ DiaSorin Liaison XL (n=93)	115.0	150.3	115.2	71.0	115.9
■ DIASource CT (n=5)	129.0	189.9	143.4	87.9	127.0
■ IDS EIA (n=7)	118.3	145.2	124.2	75.5	114.8
■ IDS RIA (n=14)	152.3	177.1	146.6	94.3	146.3
■ IDS iSYS (n=12)	134.0	153.6	133.3	84.5	129.9
■ IDS iSYS New (n=16)	100.6	125.9	98.9	63.7	102.4
■ LC-MS/MS (n=14)	111.1	136.3	113.8	71.9	107.3
■ ALTM (n=162)	117.2	148.7	117.6	73.6	116.6

\* Methods with 5 or more results returned